


Specification

Indicator Buffer Tablets for volumetric analysis

285406

| | |
|------------------------------|---|
| Physical Description: | Solid |
| Product Code: | 285406 |
| Product Name: | Indicator Buffer Tablets for volumetric analysis |
| Specifications: | Suitable as water hardness indicator: passes test Disintegration time at 25°C: ≤ 100 s |
| Hazard pictograms |  |
| WGK: | 1 |
| Storage: | RT |
| Signal Word: | Attention |
| GHS Symbols: | GHS07 |
| H Phrases: | H302 H317 H319 |
| P Phrases: | P261 P280 P301+P312 P305+P351+P338 P330 P501a |
| CS: | 38220000 |
| | |

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340