



PRODUCT CODE: 396882

Conductivity Standard 84 $\mu\text{S}/\text{cm}$ (25°C)

TARIC 3822 00 00 00

PHYSICAL DATA: liquid, Miscible with water • D 1,002 •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 2.600 mg/kg • LD L0 oral man 20 mg/kg • •

WEIGHT/VOLUME INFORMATION: 1l~1,002 kg 1kg~0,998 l

COMPOSITION:

Potassium Chloride 0,00056 mol/l

SPECIFICATIONS:

CONDUCTIVITY VALUES DEPENDING ON
TEMPERATURE:

T (°C) ($\mu\text{S}/\text{cm}$)

20,0	75,8
21,0	77,5
22,0	79,1
23,0	80,7
24,0	82,4
25,0	84,0
26,0	85,6
27,0	87,3
28,0	88,9
29,0	90,6
30,0	92,2

Uncertainty

See certificate

Conductivity at 25°C

82 - 86 $\mu\text{S}/\text{cm}$

Ed.: 3 . Vig.: 10.07.2014 .

Prod.: 396882