



X075 & X575 Egg Yolk Emulsion

Description

A pasteurised emulsion of egg yolk for use in bacteriological culture media, including Bacillus Cereus Medium PREP/MYP (LAB073), PEMBA (LAB193), Brazier's CCEY Agar (LAB160) and Perfringens Agar TSC (LAB194). It may be added directly to nutrient media for the identification of *Clostridium*, *Bacillus* and *Staphylococcus* species by their lipase and/or lecithinase activity.

Technique

Please refer to the specification sheet of the appropriate dehydrated culture medium.

Appearance

Pale yellow emulsion

Hazard Classification

Refer to appropriate SDS